

FFP2 PARTICLE FILTERING HALF MASK

CE 0598 EN 149:2001+A1:2009 FFP2 NR



This mask complies with the requirements of EU Regulation (EU) 2016/425 for Personal Protective Equipment and meets the requirements of European standard EN 149:2001+A1:2009 FFP2 NR. EU type-examination for Regulation (EU) 2016/425 by Notified Body Number: 0598.SGSFimko Oy, Takomotie 8, FI-00380 Helsinki, Finland.


USER INSTRUCTION

- ① Failure to follow all instructions and limitations could seriously reduce the effectiveness of this particle half mask and could lead to illness, injury, or death.
- ② A properly selected mask is essential. Before occupational use, the wearer must be trained by the employer in the correct use of the mask in accordance with applicable safety and health standards.
- ③ This particle half mask does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life.
- ④ Discard the mask and replace with a new one if:
 - Excessive clogging of the particle half mask causes breathing difficulty.
 - The particle half mask becomes damaged.
- ⑤ Leave the contaminated area if dizziness, irritation, or other distress occurs.

THE MAIN RAW MATERIALS

- Melt-blown non-woven fabric (PP)
- Spunbond non-woven fabric (PP)

FITTING INSTRUCTION



- ① Hold the mask in hand with the nose clip up.
- ② Position the mask under the chin covering mouth and nose.
- ③ Pull the ear loops over the ear, use the adjustable hook to fix the ear loops and adjust the position to feel as comfortable as possible.
- ④ Press soft nose piece to conform snugly around the nose.
- ⑤ To check fit, cup both hands over the mask and exhale vigorously. If air flows around nose, tighten the nose clip. If air leaks around the edge, reposition the mask/ear loops for better fit. Re-check the seal and repeat the procedure until the mask is sealed properly.

WARNING

- ① Failure to follow all instructions and limitations on the use of this product, or failure to achieve proper fit, may result in damage to your health or death or will not provide the expected level of protection.
- ② This product does not supply oxygen. Use only in adequately ventilated areas containing sufficient oxygen to support life. Do not use this respirator when the oxygen concentration is less than 19.5%.
- ③ Do not use when concentrations of contaminants are immediately dangerous to health or life. Do not use this product in an explosive atmosphere.

- ④ Facial hair, beards and certain facial characteristics may reduce the effectiveness of this respirator.
- ⑤ "NR" means non reusable, this particle filtering half mask shall not be used for more than one shift. The length of time this respirator can be used depends on contaminants present but should not exceed one shift. The respirator should be replaced sooner if breathing becomes difficult.
- ⑥ Respirator should be discarded after use or when it becomes damaged or deformed in any way; no longer forms an effective seal to the face; becomes wet or visibly dirty; breathing becomes difficult.
- ⑦ The respirators must be stored and transported in their original package and protected by the storage temperature and humidity as suggested by the manufacturer.

USE LIMITATION

- ① Do not use the mask or enter or stay in a contaminated area under the following circumstance:
 - Atmosphere contains less than 19.5% oxygen.
 - If you smell or taste contaminant.
 - For protection against gases or vapors.
 - Contaminants or their concentrations are unknown or immediately dangerous to life or health.
 - In explosive atmospheres.
- ② Do not modify or misuse the mask.
- ③ Do not use the particle half mask with facial hair or any other conditions that may prevent a good face-seal as the requirements for leakage will be achieved.
- ④ Particle filtering half mask needs to be inspected prior to each use to assure there are no holes in the breathing zone other than punctures around and staples and no damage has occurred. Enlarged holes resulting from ripped or torn filter material around staple punctures are considered damage.
- ⑤ This mask helps protect against certain particulate contaminants but does not eliminate exposure to risk or contracting disease or infection. Misuse may result in sickness or death.
- ⑥ This particle half mask marked "NR", shall not be used for more than one shift.
- ⑦ This particle half mask is not suitable for users who have damaged ears.

PRODUCT TYPE

C-Shaped
165mm×105mm C-Shaped
SNN70369C FFP2 NR



【Manufacturer】 : Shandong Shengquan New Materials Co., Ltd.

【Address】 : Diaozhen Industrial Development Zone, Zhangqiu District, Jinan, Shandong, China

【Service Hotline】 : 0086-400-777-8118

【Website】 : www.shengquan.com

【E-mail】 : sq.ppe@shengquan.com

CERTIFICATE



PRODUCT STRUCTURE

Filtration+ Bacteriostasis Dual Protection

≥94% filtration efficiency

- ① Molding Layer Non-woven fabrics
- ② ③ Filtration Layer-PP Melt-blown fabric (Double Layer Pressing)
- ④ Bacteriostasis Layer-Biomass Graphene PP Non-woven fabrics
- ⑤ Skin-friendly layer-PP Non-woven fabrics

