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Minimum door height:	
Uncrated System	2 290 x 1260 mm
Crated System	2 500 x 1500 mm
Weight distribution	700 kg/m ²
Feet supporting sub frame unit - two kinds of feet	5 x - 50 x 80 mm 4 x - 50 x 110 mm
Four round feet supporting area of base frame	100 mm in diameter
Maximum point loading	50 N/cm ²

	Weight (kg)	Center of Gravity (vertical)
Basic System	1 600 kg	997 mm from the floor
Basic System + Gatan Image filter (GIF)	1 700 kg	949 mm from the floor

Note The weight of microscope when fully installed.

4.3 Unloading Area

Make sure the unloading area is easily accessible and ready to accommodate the delivery. Keep in mind that some parts of the delivery have exceptional dimensions. Also make sure the floor can support heavy loads. For detailed information about the weights and dimension of delivered crates see Delivery Properties on page 12.

In case the delivery cannot be moved to storage or final System location straight away, ensure a dry and secure environment is available.

4.4 Transportation Route

Make sure the full transportation route from the unloading area to the System space is suitable for transport of each of the microscope parts.


The parts can be transported either in the crates or when the transport route dimensions are not sufficient the parts can be taken out of the crates in the unloading area.

Especially check all surfaces, corners and lifts (weight limit and dimensions) on the transport route.

During the transportation the System is positioned on four swivel wheels. Due to the weight of the System the floor at the Customer's site must comply to the following specification or measures need to be taken to protect the Customer floor from damaging and to transform the System safely:


- Floor should be able to support 100 N/cm²
- Floor should have a hard surface (no carpet or other soft materials)
- Transporting route should have no door steps
- Transport route should have a flat surface, no slopes are allowed

Physically moving the System and it's parts by hand is allowed, without using any vehicles to pull the System on its wheels. Transport on a fork lift is allowed when using the original transport pallet the System is delivered on.

The system width without hinges which can be removed is 124 cm.
The dimension are with the reserve size x cm on both sides.

Note If there is not enough space between the System at the highest point and the ceiling, it might not be possible to successfully install or use the System.

The rest of the supporting drawings can be found in Appendix 1: Drawings and Floorplans on page 53.

5 General Requirements of a Site

To meet the ambient conditions required for the System site a standard set-up is chosen. This set-up consists of an air handling unit and a HVAC system (Heating Ventilation and Air Conditioning). In this chapter the reference data are summarized.

The following points should be taken into account when choosing a site for the system:

- Avoid sites subjected to heavy vibrations or high acoustic noise levels, e.g. locations near elevators, trains, shipping vehicles, busy roads etc.
- Avoid sites subjected to high levels of stray magnetic fields, e.g. locations near to large motors and transformers, electric railways and tram ways etc.
- The site should be isolated from the rest of the building, preferably on a concrete block.
- Type of concrete has to be Concrete C20/C25 for safe anchoring of the system to the floor.

WARNING! Insufficient quality of the concrete floor can result in unsafe anchoring of the equipment!

5.4 Vibrations, Acoustics and Pressure waves

Disturbances due to vibrations and acoustics play a significant role. Vibration levels can be limited by proper building design in combination with vibration isolated platforms.

The floor vibrations are defined by the vibration criterion (VC) curves, by Ungar and Gordon. VC levels are defined as rms speed values, integrated over one third octave bands.

The VC-curves are to be used as guidelines only, the final room behavior with regard to vibrations will be determined during a site survey performed by trained Thermo Fisher Scientific personnel using Thermo Fisher Scientific's Site Evaluation Tool.

Depending the purpose of the microscope, time constant is involved in the vibration characteristics.

VC Curves Guidelines provided by the VC-curves can be interpreted as described below. The table lists applicable VC-curves for certain direction and system.

In case that the Microscope site does not comply with the vibration specifications, different options are available to ensure undisturbed System performance. Pre-Installation information for these Vibration Isolation Systems can be found in chapter Appendix 5: Optional Configurations on page 69.

Note For detailed information about each solution and the optimal choice for a specific situation please contact your local Thermo Fisher Scientific representative.

Vibration criteria:

VC-E	<ul style="list-style-type: none">• Floor $\geq 700 \text{ kg/m}^2$, indication > 29 cm thick• Rigid support structure with small spans or heavy foundation plate• Zone VC-E separately funded, keep around by dilation free of residual building• Technical rooms, elevators, stairways and circulation areas outside dilated zone VC-E• At the location of the dilatation disconnect channels and pipes to prevent transmission of vibration. No vibration sources in zone VC-E• Technical rooms on distance, isolate vibration sources outside zone VC-E• Zone VC-E on ground floor or 1st floor
VC-F	<ul style="list-style-type: none">• Floor $\geq 1\,000 \text{ kg/m}^2$, indication > 58 cm thick, heavy foundation plate• Zone VC-F separately funded, keep around by dilation free of residual building• Technical rooms, elevators, stairways and circulation areas outside dilated zone VC-F• At the location of the dilatation disconnect channels and pipes to prevent transmission of vibration• No vibration sources in zone VC-F• Technical rooms on distance, isolate vibration sources outside zone VC-F• Zone VC-F on ground floor



To avoid pressure waves on the Microscope Site, it must have sliding doors rather than slamming doors.

CAUTION! The required floor thickness must be calculated based on System specifications.


8 Appendix 1: Drawings and Floorplans

Note All dimensions in this section are in mm.


2510 mm is the length of the metal Frame - swivel wheel in absolute minimum length for 2680 mm the wheels are in the worst position.

Note Roof hinges need to be removed onsite.


Note Hinges can be removed (1410 mm -> 1240 mm).


Note To achieve a minimum width of 1200 mm (as indicated in the picture "dimensions without hinges and without wheels") it is needed to remove the hoist wheel and chain restrictor. Please check the Install manual 107064 how to perform this action (Chapter 2.7.4. Remove the Hoist Wheel and Chain Restrictor).



Column and Feet position with respect to the Frame

System anchoring


Center of Gravity


System Modules

HT Tank:

Note

9.1 Pre-Installation Site Survey

A pre-installation Site Survey is performed prior to tool placement to verify that the ambient acoustics, floor vibration and EMI meet Thermo Fisher Scientific specifications.

The survey is performed by a Thermo Fisher Scientific FSE who has been trained for this task.

Site survey is performed using the following kit:

- Dactron - 4035 273 20971

Additional site surveys are required when any changes are made to the site, as changes may affect the environment and system.

Note	If the site fails to meet the specifications provided in this manual the "Performance Waiver due to Failure of Pre-installation Site Requirements" must be completed and signed, see Appendix 3: Performance Waiver on page 65.
-------------	--

The system comes with an Enclosure capable of reducing the influence of acoustics and room temperature.

Note	Until Talos Arctica system is available in the SE tool, Talos F200C should be used as a replacement tool.
-------------	--

To prepare logistic departments to transport and to plan unloading of the System crates on Customer's site an FSE has to take pictures of the unloading area during site survey, except where it is forbidden by the customer.

It is needed to take pictures of:


- access road to unloading area
- unloading area where truck will be parked
- measure height of loading dock (if present)

Note	When done, the FSE has to send these pictures to the local SSOC, except where it is forbidden by the customer.
-------------	---

11 Appendix 4: System Anchoring

11.1 Talos System Anchoring

Talos F200X G2 is located inside an Enclosure. The Enclosure frame must be permanently fixed to the ground by 1 bolt.

Tilting preventions during the use of the Hoist

The customer is responsible for preparing the floor with 1 hole, this hole must have a female M16 threaded insert for tilt prevention. Bolt position with respect to frame is pictured in Appendix 1: Drawings and Floorplans on page 53.

One threaded hole with a thread depth of at least 25 mm is required. See the Floorplan for the position of the hole.

The flatness of the floor must be demonstrated with a spirit level.

- The threaded inserts must be able to sustain pull-out loads of at least 62 kN/6300 kg force and shear loads of at least 28 kN/ 2900 kg force.
 - The M16 bolt is used to prevent the system from tumbling during hoisting 310 kg outside the enclosure. In this circumstance the force on the bolt is only 800 N (80 kg).
- An anchor system that is specified to meet these requirements is sold by Hilti:
 - Adhesive item number 00241382, Hilti HIT-RE 500-SD
 - Insert item number 00258018, Hilti HIS-N M16x170

Note	The anchor systems is delivered to the Customer by Thermo Fisher Scientific. The drilling tools are not delivered. See details in Items shipped to the Customer by Thermo Fisher Scientific on page 9.
-------------	---

Note	If an S2 compliance kit is used, the Enclosure frame must be permanently fixed to the ground by 5 additional bolts. See Appendix 5: Optional Configurations on page 69.
-------------	--

11.1.1 Prepare and install Adhesive Anchors

Note	<p>Only the Epoxy Adhesive delivered by Thermo Fisher Scientific can be applied, Epoxy Adhesives from other suppliers or other Epoxy Adhesives with a different product ID from Hilti may not be applied. See Items shipped to the Customer by Thermo Fisher Scientific on page 9.</p> <p>If the Epoxy Adhesive Anchor is not delivered to the Site, contact your local Thermo Fisher Scientific representative.</p>
-------------	--

Technical parameters of the Epoxy Adhesive Anchor:

Parameter	Value
Manufacturer	Hilti
Product ID (Manufacturer)	2123403
Anchor Type	Chemical Fastener
Material Composition	Epoxy Resin
Basic Materials	Concrete (uncracked)
Temperature Range (for usage)	-40 °C - +70 °C
Storage Temparture Range	+5 - +25 °C
Thermo Fisher Scientific order code	1032987

1. The anchor holes need to be free of electrical earth. Possible earth contact can occur when the anchor makes contact with the rebar or an additional earth grid in the concrete (rebar is the steel in reinforced concrete). This unwanted contact could cause an earth loop which will be difficult to trace once the enclosure is installed. Therefore we need to check the drilled holes for earth contact after they are made. The easiest way this can be checked is with the help of a multimeter: measure the inside of the hole with respect to the earth contact to which the microscope will be connected. After the anchors are inserted in the holes it is advised to measure with respect to earth again.
2. Prepare the two component mastic gun with a new cartridge.
3. Use an M10 bolt as a helping tool to insert (and retract) the anchor.

4. Insert enough mixture in the hole to get the anchor completely surrounded by glue. (approx. 2x full squeeze).
5. Insert the anchor in the following way:
 - a. Flush to the floor or max. 5 mm below
 - b. Perpendicular to the floor
 - c. Centered in the hole
6. Remove the superfluous mixture.
7. Allow the mixture to harden for at least 24 hours.

Note **Wait at least 24 hours before the anchor is ready to be used.**

12.6 Vibration Isolation Systems

In case that the Microscope Site does not comply to the Vibration specifications (Vibrations, Acoustics and Pressure waves on page 26), different options available to ensure undisturbed System performance.


Note	For detailed information about each solution and the optimal choice for a specific situation please contact your local Thermo Fisher Scientific representative.
-------------	--

12.6.1 Integrated Vibration Isolation System (iVIS)

iVIS is a fully integrated system that reduces the sensitivity of Thermo Scientific TEM Microscopes to floor vibrations by active damping of disturbances caused by floor vibrations.

The iVIS solution will mainly damp the sensitive frequencies of the Microscope in 1-5 Hz range. The vast majority of Customer sites will meet the floor vibration specifications if iVIS is on the configuration.

iVIS is developed to actively damp the resonance modes of the column on the air isolators. It is mounted on the Microscope base in the Thermo Fisher Scientific factory and doesn't require any pre-installation actions.

12.6.2 Active Vibration Isolation solutions

Multiple solutions from different vendors are available for Active Vibration Isolation, for example anti-vibration platforms. The usual Vibration Cancellation solutions are designed to reduce vibrations over a broad frequency range.

Note	Any deviation from the specified floor flatness will result in inefficient floor vibration cancelation.
-------------	--

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